

# Appendix

## File Formats

last update: June 2, 2009 (latest changes highlighted)

The following section describes the format of the proprietary file formats used in U-view.

### Still Image File

**File Type:** binary  
**File Extension:** dat  
**Open with:** U-view

#### File Contents:

1. File header
2. if file version  $\geq 7$ : block of 'recipe' data (128 Bytes)  
if fileheader.attachedRecipeSize  $> 0$   
note: the block is of fixed size (128 Bytes!)
3. Imageheader
4. if image version  $\geq 5$ : block of markup data (128 Bytes)  
if imageheader.attachedMarkupSize  $> 0$   
note: the block is of fixed size (128 Bytes!)
5. Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

#### Description of data structures used:

##### File header

current FILEVERSION = 7

#### File header:

size in Bytes: 104

```
struct UKFileHeader{ // sizeof(UKFileHeader):104
    char id[20];           // 20 Bytes
    short size;            // 2 Bytes
    short version;         // 2 Bytes
    short BitsPerPixel;    // 2 Bytes
                           // (=16 for Sensicam, this refers to the storage,
                           // not the acquisition, which is 12 bits)
    // 6 Bytes inserted to get to next 8 byte boundary:
    // why?: because LONGLONG needs to start at 8 byte boundary
    LONGLONG spare;        // 8 Bytes
    short ImageWidth,ImageHeight; // 4 Bytes
    short NrImages;        // 2 Bytes
```

//following 58 Bytes:

//file version 7:  $\geq 1.4.3$  03/21/06 recipe sepearate: attachedRecipeSize

short attachedRecipeSize; // 2 Bytes

BYTE spare[56]; // 56 Bytes

//file version 6:  $\geq 1.4.0$  11/29/05 image sequence 'recipe'

```

short spareShort;           // 2 Bytes
BYTE SeqRecipe[56];         // 56 Bytes
};

```

### Sequencer recipe block

file version >=7:

file header may be followed by 128 Bytes of sequencer 'recipe' block  
 attachedRecipeSize=0: no attached block  
 >0 attached 128Byte block  
*attachedRecipeSize* contains count of Bytes used for data within that block

sequencer commands	arg1	arg2	arg3	
<0 do nothing				(<0 in first node: list is empty)
0 acquire	image			
1 wait	msec			
2 subtract	image=	image	-	image
3 align	dest	reference		image
4 setspin	1=spinup 0=spindown			
5 normalize difference/sum				
6 set LEEM supply absolute				
7 set LEEM supply relative				
8 LEEM preset				
9 add and accumulate	image=	image	+	image
10 accumulate arg2 into internal buffer and display that buffer / nrcycles int arg 1				

the recipe contains *further data* – in case a more detailed explanation is needed please contact Elmitec.

### Imageheader

current IMAGEVERSION = 5

4: >= 1.1.m: LEEMdata[256]

5: >= 1.5.3 03/19/08 block of markup for image tools (cross section)

**Imageheader** (*fileheader.version* >=5):

**size in Bytes: 288**

```

struct UKImageHeader{
    short size;           // 2 Bytes
    short version;        // 2 Bytes
    short ColorScaleLow,ColorScaleHigh; // from 1.3.10
                                // before 1.3.10 4 filler Bytes
    LONGLONG imagetime;   // 8 Bytes
    short MaskXShift,MaskYShift; // from 1.3.10
                                // before 1.3.10 4 filler Bytes
    //following 4 Bytes:
    //>=1.5.3 image version 5
    BYTE useMask           // >= 1.6.5
    BYTE spare
    short attachedMarkupSize //dav never contain markups
    //<1.6.5
    BYTE spare[4];         // 4 Bytes

    short spin;            // 2 Bytes
    short LEEMdataVersion; // 2 Bytes currently version 2
    unsigned char LEEMdata[256]; // 256 Bytes Overlay data
                                // 4 Bytes filler
};

```

notes:

The '**imagetime**' member of the image header structure holds the standard Windows *FILETIME* which is explained as the following:

"The *FILETIME* structure holds an unsigned 64-bit date and time value for a file. This value represents the number of 100-nanosecond units since the beginning of January 1, 1601." MS Visual C++ provides a number of functions to convert this time.

The **LEEMdata** array containing overlay data is structured as follows:

source 1, argument 1 ... source n, argument n

**source tags:**

- 0..99: 1. LEEM2000 module #  
2. followed by name  
3. followed by 1 ASCII digit identifying the unit  
    unit codes:  
    0=none,1=V,2=mA,3=A,4=C,5=K,6=mV,7=pA,8=nA,9=uA  
4. 0 to terminate the string  
5. data: 1 float (4 Bytes)

0xff: skip

- 100: Mitutoyo micrometer readout: 2 floats (x, y coordinate)
- 101: (before 1.3.10) FOV (string max 16 char's + 0)
- 102: (before 1.3.10) varian controller #1 gauge #1 value (float)
- 103: (before 1.3.10) varian controller #1 gauge #2 value (float)
  
- 104: camera exposure in ms (float)
- 105: title (string max 16 char's + 0)

new >= 1.3.10:

- 106: varian controller #1 gauge #1 label, units value (string max 16 char's+0,+string max 4 char's+0, float) 27
- 107: varian controller #1 gauge #2 label, units value (string max 16 char's+0,+string max 4 char's+0, float) 27
- 108: varian controller #2 gauge #1 label, units value (string max 16 char's+0,+string max 4 char's+0, float) 27
- 109: varian controller #2 gauge #2 label, units value (string max 16 char's+0,+string max 4 char's+0, float) 27
- 110: FOV , camera to FOV cal. factor (string max 16 char's+0+float)
- 111: phi, theta (float,float)
- 112: spin //>=1.4.1

**obsolete:**

**Imageheader** (fileheader.version <5):

**size in Bytes: 48**

```
struct UKImageHeader{
    short size;           // 2 Bytes
    short version;        // 2 Bytes
    // 4 Bytes filler
    LONGLONG imagetime;   // 8 Bytes
    long LEEMdata1_source; // 4 Bytes
    float LEEMdata1_data; // 4 Bytes
```

```

short spin;                // 2 Bytes
short spareShort;         // 2 Bytes
float LEEMdata2_data;     // 4 Bytes
BYTE spare[16];           // 16 Bytes
};

```

### **Image markup block**

image version  $\geq 7$

image header may be followed by 128 Bytes of image markup block

attachedMarkupSize =0: no attached block

>0 attached 128Byte block

*attachedMarkupSize* contains count of Bytes used for data  
within that block

This block contains info about lines and markers (letters & numbers) which the user has placed on the image.

note: the word 'markup' does not imply a similarity to *html*.

## Video File

*File Type:* binary  
*File Extension:* dav  
*Open with:* U-view  
*File Contents:* Concatenated still images.

File contents:

**Fileheader**  
if file version  $\geq 7$ : block of 'recipe' data (128Bytes)

**Imageheader<sub>1</sub>**  
**ImageData<sub>1</sub>:**  
Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

...  
**Imageheader<sub>n</sub>**  
**ImageData<sub>1n</sub>**  
Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

## Intensity data File

**File Type:** text

**File Extension:** ivs

**Open with:** U-view or simple text editor : Microsoft Word-Pad  
Notepad does not display the text correctly, Word does not save correctly unless you follow the note below:

**Note:** don't re-save the file with a formatting editor like MS Word. If this is done, the file can not be read back into U-view. You could use MS Word but must save as plain text file.

### **File Contents:**

```
UK SOFT
software      FileVersion
IRectangle    left top right bottom
StartChannel  StartChannel
DataSection   #Channels n
Time1        intensity1
....
Timen        intensityn
last_entry
```

*number of data pairs to follow  
data pairs: time, intensity  
Exponential format*

### **example:**

```
UK SOFT
software      1
IRectangle    254 174 274 194
StartChannel  0
DataSection   4
5.050000e+003 1.251472e+006
5.220000e+003 1.252496e+006
5.270000e+003 1.253216e+006
5.380000e+003 1.254112e+006
last_entry
```